

COLARIS AP[®] Technical Specifications
Myriad Genetic Laboratories, Inc. Updated: February 5, 2013

TEST RESULTS SHOULD BE USED ONLY AFTER REVIEW OF THE FOLLOWING SPECIFICATIONS:

Description of Analysis

COLARIS AP[®]: COLARIS AP[®] testing can be performed using full sequence analysis of the *APC* and *MYH* genes, together with rearrangement testing of *APC* and *MYH* by microarray-CGH analysis. Other combinations of testing can be performed as ordered on the test request form.

Full sequence analysis: Full sequence determination of *APC* is performed in both forward and reverse directions of approximately 8532 base pairs comprising 15 exons and approximately 420 adjacent non-coding intronic base pairs. Full sequence determination of *MYH* is performed in both forward and reverse directions of approximately 1608 base pairs comprising 16 exons and approximately 450 adjacent non-coding intronic base pairs. The non-coding intronic regions of *APC* and *MYH* that are analyzed do not extend more than 20 base pairs proximal to the 5' end and 10 base pairs distal to the 3' end of each exon.

Large rearrangement analysis: All exons and promoter regions of *APC* and all exons of *MYH* are examined for evidence of deletions and duplications using microarray comparative genomic hybridization analysis (microarray-CGH).

Single Site COLARIS AP[®]: DNA sequencing analysis is performed for a targeted gene region containing the specified variant in *APC* or *MYH*. Microarray-CGH analysis is performed for all requests for single site mutation analysis of a large rearrangement in *APC* or *MYH*. In some cases, long-range PCR analysis and/or sequencing of the resulting PCR product is used to detect specific, previously reported insertions in *APC* or *MYH*.

Description of Method:

Patient samples are assigned a unique bar-code for robotic specimen tracking. DNA is extracted and purified from peripheral blood samples or buccal mouthwash samples, submitted for molecular testing.

Full sequence analysis: Aliquots of patient DNA are each subjected to polymerase chain reaction (PCR) amplification reactions. The amplified products are each directly sequenced in forward and reverse directions using fluorescent dye-labeled sequencing primers. Chromatographic tracings of each amplicon are analyzed by a proprietary computer-based review followed by visual inspection and confirmation. Genetic variants are detected by comparison with a consensus wild-type sequence constructed for each gene. All potential clinically significant variants are independently confirmed by repeated PCR amplification of the indicated gene region(s) and sequence determination as above.

Large Rearrangement Analysis: Genomic DNA from patients is analyzed by microarray-CGH analysis to determine copy number abnormalities indicative of deletion or duplication mutations across the *APC* and *MYH* genes. Approximately 800 probes have been designed to interrogate all coding exons, limited flanking intron regions, and the respective promoters of *APC*. Approximately 220 probes have been designed to interrogate all coding exons, and limited flanking intron regions of *MYH*. Each probe is analyzed using

proprietary software that compares the ratio of bound patient DNA to that of a reference DNA to indicate regions of altered copy number. The microarray design includes probes to detect deletions and duplications in multiple genes tested by MGL; however, a data masking feature is used to limit the analysis only to specific genes for which testing has been requested.

Patient samples positive for deletions or duplications are confirmed by repeat microarray analysis of the genes.

Performance Characteristics:

Analytical specificity: The incidence of a false report of a genetic variant or mutation resulting from technical error is considered negligible because of independent confirmation of all genetic variants (see above). The incidence of a false report of a genetic variant or mutation resulting from errors in specimen handling and tracking is estimated from validation studies to be less than one percent (<1%). For *APC* large rearrangement analysis no false positive results were obtained through the large rearrangement testing process that uses microarray-CGH on a set of 307 individual samples. These included 7 samples that were previously examined for deletions and duplications in *APC* by Southern blot analysis and MLPA. For *MYH* large rearrangement analysis, no false positive results were obtained through the large rearrangement testing process that uses microarray-CGH on a set of 309 individual DNA samples.

Analytical sensitivity: Failure to detect a genetic variant or mutation in the analyzed DNA regions may result from errors in specimen handling and tracking, amplification and sequencing reactions, or computer-assisted analysis and data review. The rate of such errors is estimated from validation studies to be less than one percent (<1%). The analytical sensitivity of DNA sequencing performed in both directions is estimated to be >99%. The large rearrangement testing process, using microarray-CGH correctly identified all 7 samples that were positive for an *APC* large rearrangement, previously detected by Southern blot analysis and MLPA. The large rearrangement testing process, using microarray-CGH correctly identified a homozygous *MYH* deletion (previously detected by long range PCR analysis and sequencing), tested in duplicate, and a synthetic positive sample tested in four replicates, among 309 samples tested in our validation study. The synthetic positive sample was created from genomic DNA that was digested with specific restriction enzymes.

Limitations of method: There may be limited portions of *APC* and *MYH* for which sequence determination can be performed only in the forward or reverse direction. Unequal allele amplification may result from rare polymorphisms under primer sites.

MYH mutation panel analysis includes Y165C and G382D, and does not rule out the possibility of other mutations. Unequal allele amplification may result from rare polymorphisms under primer sites.

Additional information about performance characteristics of the *MYH* analysis is available in the *MYH* Technical Specifications.

TEST RESULTS SHOULD BE USED ONLY AFTER REVIEW OF THE FOLLOWING SPECIFICATIONS:

The large rearrangement analyses described above will detect deletion and duplication rearrangements involving the promoter and coding exons of *APC* and *MYH*.

These assays will not detect some types of errors in RNA transcript processing, regulatory mutations, or balanced rearrangements (i.e. inversions). Insertions that do not result in duplications will generally not be detected by microarray-CGH.

Description of Nomenclature:

All mutations and genetic variants are named according to the convention of Beaudet and Tsui. (Beaudet AL, Tsui LC. A suggested nomenclature for designating mutations. *Hum Mut* 1993; 2:245-248). Nucleotide numbering starts at the first translated base of *APC* and *MYH*.

Interpretive Criteria:

The classification and interpretation of all variants identified in the assay reflects the current state of scientific understanding at the time the report is issued. In some instances, the classification and interpretation of variants may change as scientific information becomes available.

“Positive for a deleterious mutation”: Includes clinically significant nonsense and frameshift mutations that prematurely truncate the protein.

In addition, specific missense mutations and non-coding intervening sequence (IVS) mutations are recognized as deleterious on the basis of data derived from linkage analysis of high risk families, functional assays, statistical analysis, biochemical evidence and/or demonstration of abnormal mRNA transcript processing.

Deletions and duplications of an entire exon(s) identified by Microarray-CGH analysis may also be interpreted to be deleterious. Deleterious large genomic rearrangements include single exon and multi exonic deletions and duplications that are out of frame. In frame deletions/duplications are interpreted on an individual basis and the specific evidence supporting the classification of these mutations is included in the individual patient report.

“Genetic variant, suspected deleterious”: Includes genetic variants for which the available evidence indicates a likelihood, but not proof, that the mutation is deleterious. The specific evidence supporting such an interpretation will be summarized for individual variants on each such report.

“Genetic variant, favor polymorphism”: Includes genetic variants for which available evidence indicates that the variant is highly unlikely to contribute substantially to cancer risk. The specific evidence supporting such an interpretation will be summarized for individual variants on each such report.

“Genetic variant of uncertain significance”: Includes missense variants and variants that occur in analyzed intronic regions whose clinical significance has not yet been determined. These include nonsense and frameshift mutations that occur very close to the normal stop codon.

A genetic variant of uncertain significance in *APC* is considered less likely to be deleterious if it has been observed in one or more individuals with a known deleterious mutation in the same gene.

“No deleterious mutation detected”: Includes genetic variants for which published data demonstrate absence of substantial clinical significance. Also includes mutations in the protein-coding region that neither alter the amino acid sequence nor are predicted to significantly affect exon splicing, and base pair alterations in non-coding portions of the gene that have been demonstrated to have no deleterious effect on the length or stability of the mRNA transcript. Data on polymorphic variants are available upon request.

There may be uncommon genetic abnormalities in *APC* and *MYH* that will not be detected by Colaris APSM (see **Limitations of method**, above). This analysis, however, is believed to rule out the majority of abnormalities in *APC* and *MYH*, which are believed to be responsible for most cases of Familial Adenomatous Polyposis (FAP)/attenuated FAP (AFAP) and *MYH*-associated polyposis (MAP). Data on polymorphic variants are available upon request.

“Specific variant/mutation not identified”: Indicates that specific and designated mutations or variants are not present in the individual being tested. If one (or rarely two) specific deleterious mutations have been identified in a family member, a negative analysis for the specific mutation(s) indicates that the tested individual is likely at the general population risk of developing FAP- / MAP-associated cancers.

“Positive for two MYH mutations”: Includes observations of two *MYH* mutations, or observations of two alleles of one mutation. The presence of two *MYH* mutations has been documented in recent literature to be associated with colorectal polyposis and cancer.

Causal mutations include nonsense and frameshift mutations, as well as specific missense mutations and non-coding intervening sequence (IVS) mutations recognized as deleterious on the basis of data derived from functional assays, biochemical evidence, demonstration of abnormal mRNA transcript processing and/or segregation analysis in families.

Deletions and duplications of an entire exon(s) identified by microarray-CGH may also be interpreted to be causal mutations. Deleterious large genomic rearrangements include single exon and multi exonic deletions and duplications that are out of frame. In frame deletions/duplications are interpreted on an individual basis and the specific evidence supporting the classification of these mutations is included in the individual patient report.

“Positive for two MYH mutations, clinical significance uncertain”: Includes observations of two *MYH* mutations but it cannot be determined from this analysis alone whether these two mutations are on opposite alleles. Testing one of this patient’s parents or children will determine if these mutations are on opposite alleles, thereby indicating if this patient is at increased risk of colorectal polyposis and cancer. If these mutations are on the same allele, it is currently unknown if this patient is at some measure of increased risk for colorectal polyposis and cancer.

“One MYH mutation detected, colorectal polyposis and cancer risk unknown”: Includes observations of one allele of a causal mutation. It is currently unknown whether individuals who carry a single *MYH* mutation are at some

COLARIS AP[®] Technical Specifications
Myriad Genetic Laboratories, Inc. Updated: February 5, 2013

TEST RESULTS SHOULD BE USED ONLY AFTER REVIEW OF THE FOLLOWING SPECIFICATIONS:

measure of increased risk for colorectal polyposis and cancer. Patients with one *MYH* mutation detected through the *MYH* mutation panel will automatically receive full sequence analysis of the *MYH* gene.

Change of mutation/variant classification and issuance of amended reports: Whenever there is a change in the classification of a mutation/variant within a patient's test result, an amended report will - be provided by Myriad Genetic Laboratories.